

Dry Eye Syndrome has many causes.  
And one simple solution.

For more information about Surgical Specialties™ Ophthalmic Products or to place an order, contact our customer service team at: **877.991.1110** (USA and Canada only), **+1.610.404.3520** (International), or visit [SurgicalSpecialties.com](http://SurgicalSpecialties.com)

***UltraPlug***<sup>TM</sup>  
***Punctal Plugs***


# UltraPlug™ Punctal Plugs

## Extend Your Treatment Options

**EXTENDED WEAR  
SYNTHETIC ABSORBABLE  
(PCL)**  
Evaluate within 2 months

**COLLAGEN**  
Evaluate within 1 week


**SILICONE**  
Permanent, easily removable



## UltraPlug™ Extended Wear Synthetic Absorbable Punctal Plugs

Inserted into the punctal opening to block tear drainage through the canaliculus, synthetic absorbable plugs provide extended temporary relief of dry eye symptoms. **2 plugs per package /10 packages per box.**

- Indicated for use in the treatment of Dry Eye Syndrome and the dry eye components of various ocular surface diseases
- May be used after ocular surgery to prevent dry eye related complications
- Ideal for Post-Op dry eye management of LASIK patients
- Treatment option for seasonal allergy patients experiencing dry eye symptoms
- May be beneficial for contact lens patients experiencing dry eye related symptoms
- Manufactured from E-Caprolactone-L-Lactide copolymer (PCL)
- Essentially absorbed within 2-6 months (60-180 days)




Item Number	Description
<b>EPLUG2</b>	0.2mm diameter /2.0mm length
<b>EPLUG3</b>	0.3mm diameter /2.0mm length
<b>EPLUG4</b>	0.4mm diameter /2.0mm length

## UltraPlug™ Collagen Absorbable Punctal Plugs\*

Inserted into the punctal opening to block tear drainage through the canaliculus, collagen plugs provide temporary relief of dry eye symptoms. **6 plugs per package /10 packages per box.**

- Manufactured from precisely measured premium collagen
- Used as a diagnostic assessment for permanent punctal occlusion
- Each individually packaged card contains 6 plugs



Item Number	Description
<b>C-PLUG2</b>	0.2mm diameter /2.0mm length
<b>C-PLUG3</b>	0.3mm diameter /2.0mm length
<b>C-PLUG4</b>	0.4mm diameter /2.0mm length
<b>C-PLUG-V</b>	4 Cards each of 3 different sizes - 0.2, 0.3, 0.4mm diameter /2.0mm length


\*UltraPlug™ Collagen Absorbable Punctal Plugs are not CE marked.

## UltraPlug™ Punctal Plugs

### UltraPlug™ Silicone Punctal Plugs

Inserted into the punctal opening to block tear drainage through the canaliculus, silicone plugs provide long-term relief of dry eye symptoms, and enable the eye to retain more of its own natural tears. **2 per box.**

- Low-profile cap design virtually eliminates foreign body sensation
- Flexibility of cap design facilitates easy removal with jeweler's forceps, if necessary
- Individually packaged and preloaded on a sterile, disposable inserter / dilator
- Manufactured from soft silicone utilizing a single piece mold design for greater patient comfort



Item Number	Description
<b>PPLUG4</b>	0.4mm plug / inserter
<b>PPLUG5</b>	0.5mm plug / inserter
<b>PPLUG6</b>	0.6mm plug / inserter
<b>PPLUG7</b>	0.7mm plug / inserter
<b>PPLUG8</b>	0.8mm plug / inserter

## UltraPlug™ Punctal Plugs

### UltraPlug™ Punctal Plug Educational Materials

UltraPlug™ educational materials for staff and patients in a variety of media



Item Number	Description
<b>OPG-PM-038</b>	UltraPlug™ Punctal Plug Patient Pamphlet. 50 per pack.