

EXPERIENCE EXCEPTIONAL PERFORMANCE

FROM FIRST TO FINAL STITCH



SharpPoint™ PLUS Surgical Suture


- Exceptional Needle Performance
- Thoughtfully Crafted Needle Tip
- Unique Needle Body Construction

Best In Class Performance

SharpPoint™ PLUS surgical sutures: Delivering outstanding needle performance at the first pass and superior performance over multiple passes¹

SharpPoint™ PLUS Sutures: Comparable To Ethicon™ Premium Reverse Cutting Needles

SSC™ vs. Ethicon™ PS-2 Needle Penetration¹




Up to 25% less force required versus Ethicon™ PS-2 Needles



Benchmark testing of 4-0 SSC MONODERM™ Sutures compared with Ethicon™ MONOCRYL™ Sutures.

SSC™ vs. Ethicon™ P-3 Needle Penetration¹




Up to 40% less force required versus Ethicon™ P-3 Needles



Benchmark testing of 6-0 SSC Polypropylene Sutures compared with Ethicon™ PROLENE™ Sutures.

SharpPoint™ PLUS Sutures: Exceptional Performance Versus Ethicon™ Reverse Cutting Needles

Superior to Ethicon™ FS-2 Reverse Cutting Needles¹



Up to 70% less force required versus Ethicon™ FS-2 Needles



Benchmark testing of 5-0 SSC Nylon Sutures compared with Ethicon™ Ethilon™ Sutures.

Reference: 1. Data on file. Reading, PA: Surgical Specialties Corporation; 2013. Testing performed in synthetic media.

Visit SurgicalSpecialties.com And Click "Find A Product" For SKUs And Cross-Reference Information.